



APPLICATION BRIEF - 17/0805


INDUSTRY: PHARMACEUTICAL

PRODUCT: WFI (WATER FOR INJECTION) STERILE BARRIER CLEAN STEAM BALL VALVE

MEDIA: CLEAN STEAM

IMI PBM Solution:

A major pharmaceutical company was looking for a valve to replace the traditional diaphragm valve coupled with a ball valve used as a sterile barrier for a WFI loop. PBM's valve design incorporates adjustable seats (adjust-o-seal) resulting in both an upstream and downstream seal. PBM's ability to seal at the upstream seat provided the ability to provide bubble-tight isolation with zero dead-leg at the WFI header, thus maintaining a sterile barrier. When the system is sterilized, the PBM valve is opened to introduce clean steam into the process loop. A purge port in the valve body drains any resultant condensate that may accumulate during the cleaning cycle. PBM's ability to provide an ultra-sanitary process valve that seals on the upstream seat and performs flawlessly in this application resulted in saving thousands of dollars compared to traditional methods of using a combination diaphragm valve coupled with a ball valve.



1070 SANDY HILL ROAD, IRWIN, PA 15642 USA
 1-800-967-4PBM
 WWW.PBMVALVE.COM
 INFO.PBMVALVE@IMI-CRITICAL.COM

