


Industry: Pharmaceutical
Product: Fabflex® manifold with automated sanitary ball valves and manual diaphragm valve
Media: CIP (clean-in-place) solution



Application Brief: A major pharmaceutical manufacturer required a sanitary valve manifold designed to fit in an pre-existing cabinet enclosure. The sanitary valve manifold, consisting of automated sanitary ball valves and a manual sanitary diaphragm valve is used to supply and re-circulate CIP solution to a granulator. PBM manufactured the sanitary fabflex® manifold assembly in 4 weeks so installation could take place during a scheduled shut-down. The manifold was manufactured and shipped in 3 sections (connected using a sanitary clamp and gasket) which allowed for easy on-site assembly.

APPL 0905-28